

AuditBase System requirements

Server specifications

Hardware requirements

Hardware component	Minimum	Recommended	Comment
Processor	3GHz	Dual core	
RAM	2Gb	3Gb	
HDD space	80Gb free	100Gb	
Network adapter	Required		
Other	Required		Back-up facilities
USB	1.1	2.0	

Software requirements

Software component	Description	Comment
Operating system	Windows Server2000 SP4 Windows 2003 Server SP2	All Microsoft updates should be installed

Client specifications

Hardware requirements

Hardware component	Minimum	Recommended	Comment
Processor	2GHz	Dual core	
RAM	1Gb	2Gb	
HDD space	6Gb free		
Network adapter	Required		
Graphics adapter	XGA capable of displaying 1024x768 pixels in 32bit colour		

Software requirements

Software component	Description	Recommended	Comment
Operating system	Windows 2000 Professional SP 4 Windows XP Professional SP2 Windows Vista Business Windows Vista Enterprise Windows Vista Ultimate All with SP1	Windows 2000 Professional SP4 Windows XP Professional SP3 Windows Vista Business Windows Vista Enterprise Windows Vista Ultimate All with SP1	All Microsoft updates should be installed