

# Real Ear Measurement Probe Flex

Datasheet



Auditdata

<b>Specifications</b>	Two-sided Real Ear Measurement Probe with a detachable holder.	
	Frequency Range	125Hz to 16KHz
	Maximum SPL	106dB ref, 126dB probe
	Sensitivity	-35 dB
	Equivalent Noise Level	SNR > 62dB
	Cable Length	180 cm
	Physical Weight	75 g
	Connector	8-pin Mini-DIN

## STANDARD AND OPTIONAL PARTS

<b>Standard</b>	90.3031	REM Probe Flex
	91.0304	REM-kit start pack contains:
	91.0303	Probe Tubes (40 pcs)
	91.0302	Probe Tube Guides (8 pcs)
<b>Accessories</b>	91.0303	Probe Tubes (40 pcs)
	91.0302	Probe Tube Guides (8 pcs)