

Entomed

What's New

Version SA 20X



Auditdata

The purpose of this document is to give an overview of all the new and modified features within Entomed.

New Transducer DD65v2

From now on the Entomed audiometer family SA20x is also available with the circumaural headset DD65v2 instead of the HDA200. The Headset HDA200 is discontinued since it is no longer manufactured. New software version of the EntoCard is revision 1.13.

Converting of an existing Entomed audiometer to a DD65v2 headset is basically possible, but has to be checked on a case-by-case basis and can only be done in house in combination with a full re-calibration. The basis for a successful upgrade is that the software version (EntoCard) is not lower than revision 1.12.