

# Auditbase

Office Management System

Version 6.5.0 released for Denmark, Norway,  
Sweden, Finland, Ireland and the UK

## What's New

Changes in version 6.5.0

---

**Auditdata**

## What's New

---

# Letters, attachments and reports

## User reports

- Crystal Reports runtime has been upgraded from XI to Crystal Reports 2020 engine.
- All reports created in Crystal Reports XI, 2013, 2016 and 2020 versions will work in Auditbase 6.5.0.
- You can edit the standard reports that come with Auditbase 6.5.0 using Crystal Reports XI.
- If you're thinking about which version of Crystal Reports to buy for report writing, you are definitely now advised to buy Crystal Reports 2020, which is also the cheapest.
- The report "FAMOUS - Follow-up and monitoring of new users of NHS hearing aids" has been included to support the Nottingham Clinical Trials Unit study.
- We have modified the way that separators in custom reports output to Excel to resolve issues with some regional defaults.

## Measurement and clinical

### Audiogram

- We have reinstated a dropdown for Unreliable/Reliable, now tied in with Unreliable Watermark - customers in Sweden were using the Reliable, positive value, as an audited workflow step which the watermark alone did not support.

### Audiometer and tympanometry interfaces

- The Primus Control panel has been updated to allow Primus version 5 to pass the correct transducer name from Primus Control Panel to Auditbase.
- The existing version of the Primus panel will work with Auditbase 6.5.0.
- The new Primus 5 panel will work with previous versions of Auditbase prior to 6.5.0 (other than transferring the correct transducer).

# NOAH

- Auditbase has been upgraded to use the latest Noah 4.14.2. There are no feature changes; with this version, it is mainly corrections to Noah data format translations.

## Security

### Security

- Secure AES 256 encryption will now be available on the database. This is likely to have a small impact, and is possible to have a large impact, on performance.
- There is now an option to use secure communication or encrypted messaging to and from the Auditbase client with new versions of our PAS link and updated versions of our HL7 solutions.
- The components used for encrypting Import/Export files or Auditbase passwords internally have been replaced with updated .Net encryption components.
- When you create a new user, an option will now be applied by default to force users to change their password on the first login.
- This option will be applied by Auditdata supporters to the SYSADM account for new customers.
- The Instructions For Use has been expanded as a core document stating what your responsibilities as customers are in Clinical Safety, Information Security, and other Cybersecurity aspects, including the IT environment.

## System

### Technical foundation

- Installation of Auditbase on Windows 11 is now supported with the migration to TD 7.4.5 runtime platform.